

The Waters logo is displayed in the top right corner of the slide. It consists of the word "Waters" in a white, sans-serif font, followed by a trademark symbol (TM). The background of the slide features a blue gradient with a network of white dots and lines, resembling a molecular or data structure.

Fundamentals on operation of Waters LC and control software

Waters China Ltd
Hong Kong Office

1

Content – Hardware part

The Waters logo is displayed in the top right corner of the slide. It consists of the word "Waters" in a white, sans-serif font, followed by a trademark symbol (TM). The background of the slide features a blue gradient with a network of white dots and lines, resembling a molecular or data structure.

- Introduction to Waters LC and detectors
 - Arc HPLC
 - H-Class (QSM)
 - I-Class
- Care and use of Waters LC modules
 - Routine operation and maintenance
 - Solvent consideration
- Sample preparation considerations
- Introduction to Waters LC column selection
 - Chemistry selection
 - Method development consideration
- Method transfer considerations

2

Content – Software part

Waters™

Waters CDS



Empower 3

- Configuration
 - Data management
 - System management
- Data Acquisition
 - Acquity Console
 - Inlet methods
 - Sample set
- Data Review
 - Auto and manual integration
 - Calibration curve
- Exporting data
- Report Generation

MassLynx 4.1/4.2

- Overview of MassLynx
- Data Acquisition
- Data Review
 - TargetLynx
- Data management with MassLynx
- **Waters Data Converter**

©2024 Waters Corporation

3

3

Introduction to Waters LC

Waters™

Arc HPLC



9500 PSI
HPLC

H-Class (QSM)



15000 PSI
UPLC

I-Class (BSM)



18000 PSI
UPLC

©2024 Waters Corporation

4

4

Arc HPLC

Waters™

- Quaternary HPLC pump system
 - 9500 PSI @ 5mL/min flow
 - Dual flow path design – max uptime and optional flexibility
- FTN style injection
 - Temp control 4 – 40°C
 - Max 50µL default injection
 - 96 × 2mL glass vials in two plates
- Column heater
 - Temp control 20 – 65°C
 - Max fit 3 columns



©2024 Waters Corporation

5

5

Acquity H-Class UPLC

Waters™

- Quaternary UPLC pump system
 - 15000 PSI max output @ 1 mL/min flow
 - 12000 PSI @ 2mL/min flow
- FTN style injection
 - Temp control 4 – 40°C
 - Max 10µL default injection
 - 96 × 2mL glass vials in two plates
- Column heater
 - Temp control 20 – 90°C
 - 150mm column with single guard column



©2024 Waters Corporation

6

6

Acquity I-Class UPLC

Waters™

- Binary UPLC pump system
 - 18000 PSI max output @ 1mL/min
 - 9000 PSI @ 2mL/min flow
 - Default reversed phase config
- FTN style injection
 - Temp control 4 – 40°C
 - Max 50µL default injection
 - 96 2mL glass vials in two plates
- Column heater
 - Temp control 20 –90°C
 - 150mm column with single guard column



©2024 Waters Corporation

7

7

Tubing Size, µL/foot and Usage

Waters™

Component Systems---Pump to Injector (250 µL per foot)		0.040" ID
Component Systems---Pump Plumbing (0.020" ID - 62 µL per foot)		0.020" ID
Most Alliance Tubing (0.009" ID - 12 µL per foot)		0.009" ID
Acquity Classic (0.005" ID - 3.5 µL per foot)		0.005" ID
H-Class System (0.004" ID - 2.5 µL per foot)		0.004" ID
I-Class System (0.003" ID - 1.3 µL per foot)		0.003" ID


©2024 Waters Corporation

8

8

System Band Spreading: Effect of Tubing ID and Length

Waters™



©2024 Waters Corporation

9

9

Differentiating between HPLC and UPLC

Waters™

HPLC

UPLC

- Pressure <9500 PSI or 655 bar
- Large diameter tubings
- Wide band spread
- Delay volume > 200µL
- Can run with high flow rates (≥ 1 mL/min)

- Pressure >15000 PSI or 1034 bar
- Small diameter tubings
- Tight band spread
- Delay volume within 100µL
- Generally low flow rates (<0.8 mL/min)

©2024 Waters Corporation



10

10

Detailed view – Solvent Manager, Quaternary

Waters™

- QSM-R on Arc HPLC
- For Arc HPLC path 1 is 675µL and path 2 can configure with optional mixer
- QSM on H-Class



©2024 Waters Corporation


11

11

Detailed view – Solvent Manager, Binary

Waters™

- BSM on I-class
- Pair of pump heads working in parallel




©2024 Waters Corporation

12

12

Detailed view – Flow Through Needle (FTN) Sample Manager

Waters™




©2024 Waters Corporation

13

13

Flow Through Needle (FTN) Sample Manager

Waters™



©2024 Waters Corporation

14

14


Waters recommended vials for SM-FTN

Waters™

- Routine use – use general clear glass or amber coloured vials
- When volume of samples are low, use reduced volume vials
 - Total Recovery vials
 - 300µL PP vials
- When dealing with troublesome analytes, consider HPS coated vials for maximum recovery



MaxPeak HPS




©2024 Waters Corporation

15



15

FTN Injection Mechanism

Waters™



Step 1 – Sample is aspirated into needle by the metering syringe



©2024 Waters Corporation

16

16

Flow Through Needle: Key Characteristics

Waters™

- Multiple needle and loop options, typically has a wider injection range
 - Arc HPLC injection volume range: 0.1-50 μ L standard configuration
 - UPLC injection volume range: 0.1-10 μ L standard configuration
 - Can go up to 1000 μ L with extensions (not recommended for UPLC)
- Needle in flow path design
 - No wasted sample; all volume sample injection
 - Interior of needle flushed with gradient for lowest carryover performance
- Wide range of plates and vials
 - FTN does not have a piercing needle, so some adhesive/foil plate caps can lead to clogging/blocking of needle
- Simplest mechanism for easiest usage

©2024 Waters Corporation


17

17

Column management – Arc HPLC

Waters™

- Forced air column heater
- Fit 3 columns
 - 1x 250mm
 - 1x 150mm, 1x 50mm
- Temp control 20 – 65°C
- Optional column selection valve for quick column switching (3 selections)
- Clip column onto holder with space between ferrule and column body
- Column connection with finger tight fitting, no tools required.



©2024 Waters Corporation

18

18

Column management – Acquity UPLC CH-A

Waters™

- Temperature range: From 5°C Above Ambient to 90°C
- 2.1 to 4.6 mm, up to 150 mm with filter or guard column
- Active pre-heating with finger tight fitting
 - Improves RT stability
 - Reduces backpressure



©2024 Waters Corporation

19

19

Column management – Acquity UPLC CM-30A

Waters™

- Column manager with possibility of using 4 x 50mm or 2x 150mm long columns
- Programmable column switching (max 4 selections)
- Temp control 4 – 90°C, two independent temperature zone
- Use Active Pre-Heaters for best UPLC peak results



©2024 Waters Corporation

20

20

Tightening Procedure for 3-Piece Active Pre-Heater

Waters™

- Loosen the Stainless Steel Cap Nut from Gold Fitting
- Finger tighten the Gold Fitting with ferrule into column detail
- Engage the column hardware, finger tighten column onto Gold Fitting
- Finger tighten the stainless steel Cap Nut onto Gold Fitting



Pictured Above: Active Pre-Heater prior to modifications made for "Plus" Instruments

©2024 Waters Corporation



21

21

Dead volume - fixed fittings

Waters™

- Improper fittings results in void or dead volumes
- Causes bad peaks, shouldering peaks, leakages
- Waters ferrule depth is not the same with other vendors
- Fixed fittings are generally not reusable cross vendors



©2024 Waters Corporation

22

22

Waters LC detectors

Waters™

Main focus

- 2998/Acquity eλ – PDA, 3D UV
- 2489 – TUV, dual 2D UV
- 2475 – FLR
- 2414 – RI
- 2424 – ELS

Discuss next time

- QDa (SQ-MS)



©2024 Waters Corporation

23

23

Waters 2998/Acquity eλ Detector

Waters™

- Both 2998 and Acquity eλ supports 3D UV data collection
 - 2998 on HPLC, Acquity eλ on UPLC
- 190 – 800 nm , 1.2 – 12 nm bandwidth
- 80Hz data collection
- New lamps warranty: 2000 hours or 1 year
- Under Empower control, supports one (1) 3D channel + Eight (8) 2D channels data collection simultaneously
- Flowcell not interchangeable between 2998 and Acquity eλ
- Analytes must have chromophore
- Best to keep responses below 1AU

©2024 Waters Corporation

24

24

2998 PDA

Waters™

Deuterium
lampCuvette
holder

Flow cell


©2024 Waters Corporation

25

25

Acquity PDA

Waters™

Deuterium
lampBack pressure
regulator

Flow cell

©2024 Waters Corporation

26

26

Waters 2489 dual channel UV detector

Waters™



Deuterium
lamp

Cuvette
holder

Flow cell

©2024 Waters Corporation


27

27

Read lamp energy

Waters™

- Good lamp + clean flow cell + good bench optics will produce value >20,000 at 230nm
- Must turn on lamp for 60 minutes before performing test
- Run 0.5mL/min Methanol through flow cell when testing
- WKB68243
- If reading is <20,000 then change lamp and try again
- If new lamp still gives low reading, clean flow cell or call for service
- WKB48267



©2024 Waters Corporation

28

28

Read lamp energy

Waters™

Console for System Acq PDA ELS on Node M01KMBH - [PDA Detector]

Control Configure Maintain Troubleshoot Help

System
Quaternary Solvent Manag
Sample Manager FTM
Column
PDA Detector
ELS Detector
Column Manager
Plots
Maintenance Counters
Logs

conditions

Shutter Closed

PDA Max Plot

AU

Read Energy

Exposure time: 5.0 msec

Shutter Position: Open Erbm

Continuous read

Subtract dark current

Energy Value

Wavelength: 185.8255 nm

Value: 225.4113 counts

Results

Energy (counts)

Wavelength (nm)

Click on Wavelength.

Click Read

Select as close as possible to 230 nm.

Results

Energy (counts)

Wavelength (nm)

X Axis: Wavelength Pixel

225.4113
231.776
232.5983
234.1408
235.3229
236.5063
237.6886
238.8729
240.0562
241.2396
242.4239
243.6082
244.7925
245.9768
247.1612
248.3465
249.5308
250.7151
251.9024
253.0877
254.274
255.4594
256.6457
257.832
259.0193

©2024 Waters Corporation


29

29

High sensitivity flow cell option

Waters™

- High sensitivity flow cell path length is 25mm vs standard analytical at 10mm
- Theoretically provides 2.5X sensitivity
- Will reduce resolution between closely eluted peak due to increased system volume and larger bandspread now
- Very expensive >HK\$40,000



©2024 Waters Corporation

30

30

Waters 2475 Fluorescence Detector

Waters™

- Excitation: 200 ~ 890nm
- Emission: 210 ~ 900 nm
- Up to four 2D channels or one 3D channel
 - 10 points/s in single wavelength mode
 - 1 point/s in 3D mode
- Illumination by Xenon lamp
 - Warrantied at 2000 hours or 1 year
- Flow cell:
 - 13µL
 - 145PSI or 10 bar. Limit to below 5mL/min, or check backpressure when using
 - Expensive, > HK\$45000

©2024 Waters Corporation

31

31

Waters 2475/Acquity FLR operation consideration

Waters™

- Waters FLR can operate in Emission or Energy mode
- WKB118420
- Emission mode
 - Is normalized to standard water reference
 - Peak response not usually affected by gain settings
 - Compensate variation from degradation of bench optics and lamp
 - Allows more comparable result with other 2475 or Acquity FLR
 - WKB193270
- Energy mode
 - Is the direct output from PMT
 - Influenced by gain settings
 - Usually used with legacy established methods
 - WKB29076

©2024 Waters Corporation

32

32

Waters 2475 FLR Detector

Waters™



Xenon lamp

Flow cell

©2024 Waters Corporation


33

33

Waters 2414 Refractive Index Detector

Waters™

- Detection is based on a change of optical refraction in sample cell compared to reference cell
- Generally useful for most UV transparent analytes
 - Sugars, alcohols, fatty acids etc
- Used only with isocratic method
 - Cannot run gradient, must use premixed mobile phase
 - Changing gradient, back pressure causes undesirable baseline drift
- Generally wider peak width, less sensitive and selective



©2024 Waters Corporation

34

34

Waters 2414 Refractive Index Detector

Waters™

- Flow cell backpressure rated at 100PSI
 - Further protected by pressure relieve valve
 - Contains reference and sample cell in one unit
 - Reference path content not change during sample acquisition
- Flow rate: 0.1 to 10.0 mL/min
- Internal oven: 30 to 55 °C
 - Oven takes a long time to reach set temperature
 - Turn on when reach lab if to be used later

©2024 Waters Corporation


35

35

Waters 2424 Evaporative Light Scattering Detector

Waters™

- Detection based on light scattering of very fine solid particles
 - Not gaseous
 - Nitrogen helps to blow dry the mobile phase
- Operation principle:
 - Nebulization
 - Desolvation / Evaporation of the mobile phase
 - Detection: diffusion chamber



Column outlet



©2024 Waters Corporation

36

36

Waters 2424 operation consideration

Waters™

- Evaporative Light Scattering Detector is almost a universal detector and is destructive
- Detection of light scattered from particles of analyte formed by nebulization and evaporation of the mobile phase
- Compounds have to be less volatile than mobile phase
- The sensitivity of the ELSD (nanogram level, low ppm level) is better than the refractive index detector, and can be used with gradient elution
- Non-linear calibration curve
- Use of LCMS compatible MP is highly recommended
- Siphon drain tube usually requires a column of liquid for ELS to perform best

©2024 Waters Corporation

37

37

Waters™

Routine Operation Considerations

38

Turning on PC and starting Empower 3

Waters™

- If all LC and computer are turned off, switch on PC first
- Login to Windows and wait for 5 minutes
 - Empower 3 will then be operational because it takes time to start Oracle database
- After logged into Windows, turn on all LC modules
- Switching on 2998/2489/2475/2414 from powered down state will ignite lamps
- All modules will initialise and perform startup check
- If your PC is running MassLynx as acquisition software, just turn on PC, switch on modules, and run MassLynx after logging into Windows

39

General LC operation – shutdown operation

Waters™

- Flush buffered mobile phase from system and detector with purified water
 - For Binary Pump, reserve either A1 or A2 for purified water
 - For Quaternary Pump, reserve A for purified water
 - Run separate method for flushing pathway and column
 - E.g. 0.5mL/min 90:10 water:methanol for 20 minutes, 50:50 water:methanol for 20 minutes, then 100% methanol for 10 minutes
 - Substitute methanol with acetonitrile is possible
 - Consult user manual for your column for best storage solvent
 - Can maintain pump at very low flow (0.02 mL/min) if will be using on next day
- Turn off column heating
- Store flow cell in methanol or acetonitrile
- Turn off if going for a long holiday
- Cap the inlet and outlet

40

Solvent manager – usage considerations

Waters™

- All mobile phase sinkers must be immersed in mobile phase at all time
- Mobile phase bottles must be covered appropriately at all time
- Prime all solvent lines and ensure no bubbles in it
 - Solvent refresh: >3 minutes
 - Solvent changeover: ≥5 minutes
- Use proper seal wash – 90% water + 10% methanol/isopropyl alcohol
 - Set frequency to 2 – 5 mins
 - This is compatible for all reversed phase methods
- Prime seal wash to remove air bubbles. Use syringe to pull seal wash solvent through when necessary
- Online degasser need time to draw vacuum, turn on Solvent manager upon arriving in lab to reduce startup wait time

©2024 Waters Corporation

41

41

Solvent manager – usage considerations

Waters™

- Respect pH range of instrument, 1 – 12.5
 - Do not use mineral acids (HCl, HF, HNO₃) and strong bases (NaOH, KOH) as mobile phase
 - H₂SO₄ can be used at very low conc., ~ 0.005N
 - Final operational pH must be compatible with your LC column
- Use Type 1 water for preparing aqueous buffers
 - Replace aqueous mobile phases daily or at least weekly
 - If water is your mobile phase, replace daily
- Use only HPLC grade organic solvents or above
 - Or ACS grade if not available
- When using salts, use HPLC, LCMS or ACS grade
 - Generally highest quality that you can acquire, of reputable source
 - Concentrated buffer (< 0.1M) can be used but generally please avoid

©2024 Waters Corporation

42

42

Solvent manager – Recommended solvents/modifiers

Waters™

- Waters HPLC/UPLC System recommended solvents
 - Methanol, Acetonitrile, Isopropanol and their water mixtures
 - THF (replace all PEEK tubings with Stainless steel, THF/Hexane modification kit, p/n 205000464)
- Sample diluents (in addition to the solvents listed above)
 - Dimethyl sulfoxide (DMSO)
 - Dimethylformamide (DMF)
- Additives / Modifiers
 - 0.2% formic acid, 0.1% trifluoroacetic acid (TFA), 0.1% triethyl amine (TEA), phosphate buffer, ammonium bicarbonate, ammonium hydroxide, ammonium acetate, Ethylenediaminetetraacetic acid (EDTA)
- Cleaners
 - Phosphoric acid (~30%), Sodium hydroxide (~1M)
 - Ask before using these

©2024 Waters Corporation

43

43

Solvent manager – usage considerations

Waters™

- Filter all buffered mobile phase through 0.2 µm membrane filter or risk clogging your LC and column
 - Filtering also degas them
- When preparing isocratic premixed aqueous organic mobile phase, degas before loading onto the LC
- Degas mobile phase before use produce flatter baseline – sonication or vacuum filtration
- Don't top up aqueous mobile phase, especially during acquisition
- Store mobile phases in clean bottle
 - For buffered mobile phase, consider amber bottles to reduce microbial growth
 - Label all mobile phase properly, including content and date prepared
 - Cover exposed mouth or openings properly

©2024 Waters Corporation

44

44

Solvent manager – filtration set

Waters™

- Using filtered mobile phase is very important
- Laziness will cost you instrument downtime and column lifetime
- Filtration will also degass your mobile phase and improves chromatographic performance



©2024 Waters Corporation

45

45

Sample manager – common issues

Waters™

- Issues are generally physical or LC method related issues
- Physical
 - Bent needle: misalignment when drawing sample, vial inserts, user ignorance
 - Blunt needle: due to thick septum or non-preslit caps
 - Lost prime in solvent tubing serving FTN: trapped air in tubings
 - Worn injection valve: high usage or corrosive solvent
 - Syringe leak: high usage or corrosive solvent
- Method related
 - Carryover due to insufficient wash time/bad solvent choice
 - Wrong solvent in purge
 - Sample too concentrated

©2024 Waters Corporation

46

46

SM-FTN – solvent consideration and tubing management

Waters™

- SM-FTN served with two solvent tubing:
 - Sample manager purge : 90% water + 10% methanol
 - Sample manager wash : 50% methanol
- Sample manager purge solvent line for auto dilution and hydraulic action of SM-FTN
 - Usually 90% aqueous
 - When using HILIC, use initial condition ratio in gradient table
- Sample manager wash for washing exterior of sampling needle
 - Usually 50% methanol, choose a solvent that will dissolve your sample/analytes readily
- Both solvent tubings must be filled and without air bubbles for reliable performance

©2024 Waters Corporation

47

47

Sample Manager – extending lifetime of needle

Waters™

- Use preslit caps
 - 186000305
- Do not use crimp style GC-caps – too thick to puncture
- Do not use vial inserts, they will damage needle
 - For low volume sample, consider 300 µL PP vials 186002639
 - Or consider Total Recovery vial 186000384C
- Do not use buffered solution for washing needle.
- Using low concentration of acid such as 0.1% FA is fine



©2024 Waters Corporation

48


48

Sample Manager – loading sample plates

Waters™

- When required to load or remove sample plates, pull the sample tray handle out to do so
- Failure to do so may result in hitting autosampler arm, bending needle or misalignment

Figure 2-14: Loading sample plate onto sample tray



©2024 Waters Corporation

49

49

Sample Manager FTN usage consideration

Waters™

- Temp control 4 – 40°C, commonly used at 10°C when temperature control is required
 - Advise not to cool for extended time, use refrigerator for sample storage
- Do not overtighten caps
- Injection more than 50 µL on Arc HPLC is possible by installing larger extension loop
- Air gap can be used when injecting very large or small volume to protect sample during injection and reduce sample dispersion
- Arc HPLC control column temperature through SM-FTN page

©2024 Waters Corporation

50

50

Sample Manager FTN usage consideration

Waters™

- Priming Recommendations →
- Selecting Wash Solvents
 - Sample should be very soluble in the wash solvent
 - Purge solvent
 - does not contact sample
 - should be similar to initial mobile phase conditions, 10% MeOH is very commonly used
 - should not contain additives

To Refresh Solvents

Prime

☒ Wash Solvent: 15 sec

☒ Purge Solvent: 5 cycles

OK Cancel

For New Solvents

Prime

☒ Wash Solvent: 120 sec

☒ Purge Solvent: 50 cycles

OK Cancel

©2024 Waters Corporation

51

51

Alternative to APH/Arc Fingertight fitting

Waters™

- APH are expensive, >HK\$8000
- Improper handling of APH results in stuck ferrule, bent and deformed tubing
- Use PEEK tubings, Arc HPLC can use 0.010 or 0.005" i.d. tubings
- H-Class and I-Class must use stainless steel tubings for >6000 PSI application
- Use either fixed ferrule tubings or other finger tight types
- Waters sell MarvelX UPLC tubing finger tight to 18,000 PSI and reusable for over 100 times



MarvelX UPLC tubing

©2024 Waters Corporation

52

52

Acquity In-line filter – protects your precious columns

Waters™

- Houses a 0.2µm frit – effectively trap any particulates in mobile phase or sample that may otherwise clog columns
- Protects column from premature failure due to overpressure
- Frits are easily replaceable and cost ≤ HK\$100
 - May need to replace nut also when frit cannot be separated from it
- Each in-line filter are configured for that brand of column, using on other brand may result in leakages or bad peak shape



©2024 Waters Corporation

53


53

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Flow cells are very expensive, >HK\$40,000
- Inject only filtered samples and use only filtered mobile phases
 - Clogged flow cell or outlet tubing will result in increased back pressure causing flow cell to break
 - Filter all mobile phase or sample with 0.2µM membrane filter
 - When using concentrated buffer MP, ensure the salt do not precipitate out when switching to high organic which may cause blockages
- Flow cell will withstand 1000 PSI back pressure

Flow Cell, ACQUITY Arc Detector (# 205001552)



©2024 Waters Corporation

54

54

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Powering on instrument using front switch will always turn lamp on
- Understand which wavelength is the cutoff of your MP or modifier
 - E.g., TFA ~210nm, Acetate ~ 240nm, Methanol ~205 nm
- Collect sufficient points within a peak for best peak shape
 - 25-50 points per peak
- Check lamp energy regularly
- Turn off lamp after completing sample sets
- Flush buffered mobile phase from PDA flow cell after use
- Store PDA flow cell in Methanol or ACN
 - For long holiday, cap inlet and outlet

©2024 Waters Corporation

55

55

Cleaning PDA/TUV flow cell

Waters™

- Dirty flow cell issues:

Lowered intensity	High baseline noise
Baseline fluctuations	Calibration failure

- General cleaning procedure
 - Connect system with union
 - Flush system and flow cell with water, 30 mins
 - Prime one or all line with magic mix (25% each: MeOH, IPA, ACN, Water w/ 0.1% FA)
 - Set column temperature to 40 °C and route flow path through column heater/manager
 - Set flow to 0.2 mL/min, 60 mins
 - Flush system and flow cell with water until neutral, set flow 0.2 ~ 0.5 mL/min or do not exceed 1000 PSI
 - Setup system and try run
 - If blockage of flow cell suspected, try backflushing
 - WKB48267

©2024 Waters Corporation


56

56

Waters 2424 operation consideration

Waters™

- Nebuliser is generally set to cooling temperature to prevent introducing too many particles or analytes into the detection chamber
- After completing samples, run 50% methanol to clean the nebulizer
- Blow dry the ELS detector for another 1 hour with nitrogen
- If blockage in nebulizer is suspected, remove and sonicate with 50% methanol



©2024 Waters Corporation


57

57


Waters 2424 – removing nebuliser

Waters™

Quick-disconnect
tubing fitting at
twelve o'clock position



Removing the nebulizer:



Pins inside the desolvation chamber

©2024 Waters Corporation

58


58

Using Empower 3 for data acquisition and data analysis

59

Empower 3 – login

- Look for shortcut on desktop or Start menu
- Empower 3 requires password to unlock
 - Username: system
 - Password: manager
- Do not change password, if you forgot it and call for repair service, Waters will charge for that
- If you rebooted the computer, wait for 5 minutes before login because Oracle database is starting in background



60

Empower 3 – main screen

Waters™

- Empower Help files
- Empower information
- Administration tasks
- Data Acquisition and instrument control
- Data review and analysis
- Exit from Empower




©2024 Waters Corporation

61

61

Empower 3 – Bringing up console

Waters™



©2024 Waters Corporation

62


62

Waters 2998/Acquity eλ Detector operation consideration


Waters™

- Always check Lamp On box for acquisition method
 - Uncheck Lamp On for shutdown method
- Collect 3D data saves your time to rerun sample at specific wavelength and can provide spectra information
 - Use desired range, or check literature
 - Use highest resolution for best quality data

Figure 1-4: Benzene spectrum at different resolutions



©2024 Waters Corporation




63

63


Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Use sampling rate to maximise SN ratio and collecting sufficient points across a peak for confident quantitation
 - Around 20 – 30 points per peak required
 - HPLC methods use 10 or 20 pts/sec
 - UPLC methods use 20 or 40, 80 when using ballistic gradients



©2024 Waters Corporation



64

64

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Use sampling rate to maximise S/N ratio and collecting sufficient points across a peak for confident quantitation
 - Around 25 – 50 points per peak required
 - HPLC methods use 10 or 20 pts/sec
 - UPLC methods use 20 or 40, 80 when using ballistic gradients
- Using sampling rate that is too high may result in more noise, decreasing S/N ratio

Figure 1-6: Example of how baseline noise increases with higher sampling rates





Figure 1. Impact of different detector settings on the peak shape for acenaphthene. The number in parenthesis is the filter time constant, in seconds.

©2024 Waters Corporation


65

65

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Filter time constant
 - Noise reduction calculation performed in instrument
 - Default is normal
 - Can turn off
 - Does not affect peak area
 - For max resolution, use Fast
 - For max sensitivity, use Normal
 - When set to Fast, baseline noise greatly reduced but will shorten or broaden peaks
 - When set to Slow, removes less baseline noise, peak shape maintained narrower, small peaks more difficult to pick out from baseline
- WKB90261 and Tip207



©2024 Waters Corporation


66

66

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- D2 lamp produced peak at 370 and 656 nm
- Enabling will tell detector to ignore signal from these diodes
- When you disable this parameter, the detector reports the signal from the photodiode at 656 nm
- Disable this parameter only if you are working with compounds that absorb in the 656-nm range
- WKB44903



©2024 Waters Corporation


67

67

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Max eight other 2D channels
- These functions are available
 - Absorbance
 - Absorbance –Compensation
 - Absorbance –MBF
 - Max Plot
 - Difference
 - Sum
 - Ratio
- Resolution of 4.8nm works for many analytes in both HPLC & UPLC
 - Optimisable further for sensitivity and S/N ratio



©2024 Waters Corporation

68


68

Waters 2998/Acquity eλ Detector operation consideration


Waters™

- If monitoring at analyte λ_{max} , higher bandwidth generally reduced peak height

Figure 1-8: Resolution comparison for anthracene



©2024 Waters Corporation




69


69

Waters 2998/Acquity eλ Detector operation consideration

Waters™



- WKB67349



©2024 Waters Corporation


70

70

Waters 2998/Acquity eλ Detector operation consideration

Waters™

- Usually will add a Lamp Off event
- The Time is your actual run time + 1 minute
 - i.e. if your acquisition time is 10 mins, set to 11 mins
- In the event of Empower disconnected from the running instrument, the 11th minute will trigger a lamp off
- This protects the lamp and detector
- But if the time was not changed when sample was required to run pass 11th minute, no meaningful data will be acquired




©2024 Waters Corporation

71

71

ACQUITY UPLC® Independent Sampling Rates and Filter Time Constants

Waters™



Benefit: High Sensitivity
with Optimal Resolution

©2024 Waters Corporation


72

72

UV Spectra Resolution (Bandwidth) Effects on Peak Identification

Waters™

- High spectral resolution leads to easier ID of components
- Lower resolutions shift UV maxima
- High Confidence in Spectral Homogeneity Checks (Peak Purity)
- High confidence in UV/Spectral Matching



Spectral Resolution
1.2 nm vs. 3.6 nm

©2024 Waters Corporation

73

73

Waters 2475/Acquity FLR principle and consideration

Waters™

- Most molecule do not fluoresces, derivatisation may be required to detect certain compounds
 - E.g. amino acid, aflatoxins etc
- Principle: in a dark environment, excite everything with some light, then measure emitted light
- Fluorescence is proportional to :
 - Source intensity i.e. more intense source = larger signal
 - Volume of flow cell i.e. longer path length means more photons captured = larger signal
- High temperature decreases fluoresce intensity
- Response is only linear at low concentration, analyte dependent
- Normalisation of FLR needs to be performed regularly so to ensure results are comparable. WKB193270

©2024 Waters Corporation

74

74

Waters 2475/Acquity FLR operation consideration

Waters™

- Waters FLR detectors flow cell can only withstand 500 PSI back pressure
- Must ensure FLR is the last in system connection
- FLR flow cells are very expensive, ~HK\$45,000
- Warmup FLR for one hour before starting data acquisition
- When Xenon lamp is operating, temperature is very high and will cause burns when mishandled
- If 2475 failed to normalise during startup, try flushing flow cell clean and then run again
- Clean water and fresh mobile phase produce cleaner baseline and better S/N ratio
- FLR PMT will get saturated at very high concentration

©2024 Waters Corporation

75

75

Waters 2475/Acquity FLR usage consideration

Waters™

- Saturation of PMT, due to:
 - Incorrect gain setting
 - Analyte concentration too high
- Use timed segment gain settings
- Dilute sample
- Use another excitation wavelength
- WKB17612




Figure 3-19 Gain Set too High (Fluorescence Pinned to -9999.9 EU)

©2024 Waters Corporation


76

76

Waters 2475/Acquity FLR operation consideration

Waters™

- When operating in 2D mode, 4 channels are available
- Only a single 3D mode is available
- 2D and 3D cannot be simultaneously used in acquisition
- Supports the following data mode in 2D:
 - Emission
 - Difference A-B
 - Difference B-A
 - Max Plot
 - Reference Energy



©2024 Waters Corporation


77

77

Waters 2475/Acquity FLR operation consideration

Waters™

- Noise filter:
 - None
 - Hamming, digital
 - RC, analog
- For general purpose, Hamming filter and 1 sec is used
- Effect of:
 - Lower time constant
 - Remove less baseline noise
 - Small peaks hard to discriminate from baseline noise
 - Narrower peaks
 - Higher time constant
 - Remove more baseline noise
 - Broader and shorter peaks



©2024 Waters Corporation


78

78

Waters 2475/Acquity FLR operation consideration

Waters™

- Under 3D mode, acquisition scan type can be configured as Excitation or Emission
- The instrument will then perform a scan in the other mode across the specified wavelength range
- PMT gain depends on your analyte and concentration, adjust as required
- Noise filter hamming and 3 sec enough for general use



©2024 Waters Corporation


79

79

Waters 2475/Acquity FLR operation consideration

Waters™

- Events should be enabled for turning off lamp if computer lost connection with running instrument
 - Add an event of run time + 5 mins to turn off lamp
- Can also be used for optimizing gain across detection window
- WKB15462



©2024 Waters Corporation

80


80

Waters 2475/Acquity FLR operation consideration

Waters™

- Ample time should be allowed between changing of gain value to stabilize baseline
- Peaks that are too close cannot be used with time segment, compromise is needed

Gain optimized chromatogram:



©2024 Waters Corporation


81

81

2414 Detector - General

Waters™

- Method type can be acquisition or purge
- When selected to use purge, RI purge valve will open
- Mobile phase will flow through both reference and sample cell
- Used when conditioning the entire system
- All other functions will be locked out
- Unless mobile phase is very precious, recycle valve can be unchecked




©2024 Waters Corporation

82

82

Waters 2414 RI Detector

- When method type is in acquisition mode, all settings can be altered
- Sampling rate use 5 or 10, method dependent
- Filter time constant use Normal for most application
- Polarity will affect peaks upright positioning, if they come out as negative peaks then change to negative
- Should always autozero on inject start
- Recycle off
- Set target flow cell temperature, usually 40 is fine
- External column temperature is not required to set here




©2024 Waters Corporation

83

83

Arc HPLC QSM-R - General

- 1. Solvent information only, no actual physical performance affected
- 2. Set pressure limit according to your column
 - HPLC column: generally ~6000 PSI
 - UHPLC column: generally up to 9000 PSI
 - Other specialty column: Check user manual
- 3. Gradient timetable, set according to desired gradient
 - Gradient curve available for method transfer or delicate control of gradient profile
 - Comment for information
 - Do not adjust values in Advanced tab



©2024 Waters Corporation


84

84

Arc HPLC QSM-R – General/Auto Blend

Waters™

- Arc HPLC QSM can be programmed to perform automatic mixing of mobile phase to produce different concentration of salt and pH
- Useful for unattended method development or column screening
- Refer to 720006557 application note




©2024 Waters Corporation

85

85

Arc HPLC QSM-R – Gradient curve

- Generally useful for delicate method development
- When used, transferring method to another LC may be challenging
- 6 is always linear ramp
- 1 is immediate ramp at beginning
- 11 is immediate ramp at end
- WKB3489




©2024 Waters Corporation

86

86

Arc HPLC QSM-R – Gradient SmartStart

- Do not adjust accelerate speed
- Gradient starts is beneficial for method transfer
 - At injection: system inject samples and starts gradient
 - Before injection: system starts gradient and then injects, peaks eluted earlier
 - After injection: system gradient pauses after injection, then begin gradient
- Does not adjust gradient table, less impact on regulatory standpoint




©2024 Waters Corporation

87

87

Arc HPLC QSM-R – Data

- Keep System Pressure data is always recommended
 - For diagnostic use
 - For keeping track of column backpressure
- Keeping too much data is not recommended, will cause too many channels to be recorded



©2024 Waters Corporation


88

88

Arc HPLC FTN-R – General

Waters™

- Wash Solvent
 - Choose something that will dissolve all your analytes
 - 50% methanol is a good choice
 - Can choose to wash needle before it sit into seal, or after sitting and inject the samples
 - Wash longer for sticky samples
- Purge solvent fills the syringe, part of hydraulics in injection system
- Selecting whichever solvent will not impact system performance, only information only



©2024 Waters Corporation


89

89

Arc HPLC FTN-R – General

Waters™

- Active preheater – generally not available in Arc HPLC config
- Loop ahead and loop offline used for high throughput, don't use for ordinary LC analysis
- Column position for choosing Column if system configured with column selection valve
- Sample compartment temperature can be controlled between 4 and 40 °C
- Column temperature configurable depending on column compartment
 - CH/C 20 – 65°C
 - CH-A 20 – 90°C
 - Alarm band is optional



©2024 Waters Corporation


90

90

Arc HPLC FTN-R – General/Advanced

Waters™

- Do not adjust settings related to:
 - Syringe draw rate
 - Needle placement
 - Air gaps
 - Unless you know what you are trying to do
 - This carries a risk to contaminating entire injection system and bending needle
- Auto Addition can be used if needed to perform online mixing of samples and reagent
- No injection can be used to diagnose system or injection carryover. This option allows for recording of data



©2024 Waters Corporation


91

91

Arc HPLC FTN-R – Dilution

Waters™

- Arc HPLC FTN-R can perform online dilution
- Recommended to dilute offline



©2024 Waters Corporation


92

92

FTN – Wash Mechanisms


Waters™

Pre-Inject Wash



- Exterior of needle is washed
- Needle **NOT** in the seal, **NOT** under pressure during washing
- Wash solvent flows from bottom of port and exits through top
- Occurs before injection is started

Post-Inject Wash



- Exterior of needle is washed
- Needle is in the seal at high pressure during washing
- Wash solvent flows from bottom of port and exits through top
- Occurs while injection is running

©2024 Waters Corporation


93

93

2414 Refractive Index Detector

Waters™

- Method type:
 - Acquisition: acquire data
 - Purge: do not acquire data, will purge reference cell only. Choose this for unattended instrument conditioning method
- Sampling rate:
 - Choose between 2, 5 and 10
 - Aim between 25 to 50 points across peak
 - Usually HPLC methods use 5
 - UHPLC methods use 10



©2024 Waters Corporation


94

94

2414 Refractive Index Detector

Waters™

- Filter time constant:
 - Normal is fine
 - Other settings requires optimization
 - RI methods usually too long to make it worthwhile to optimise this setting
- Polarity:
 - Positive or negative
 - If all peaks came out as negative peaks, then switch polarity. Vice versa
- Autozero:
 - Should be checked
- Divert valve to Recycle:
 - Generally not used to rid the possibility of contaminating mobile phase reservoir



©2024 Waters Corporation


95

95

2414 Refractive Index Detector

Waters™

- Flow cell temperature:
 - Generally do not affect peak shape
 - Use 40°C or 5 °C above ambient
- External column
 - Not used
- Comment
 - Type in notes for reminder, will not affect acquisition




©2024 Waters Corporation

96

96

2414 Refractive Index Detector

- Data channel
 - Must be checked, else no data points recorded




©2024 Waters Corporation

97

97

2414 Refractive Index Detector

- Analog out is not usually configured
- Can be safely ignored for most of the time



©2024 Waters Corporation


98

98

2414 Refractive Index Detector

Waters™

- Events control is beneficial when running for extended time and unattended shutdown of instrument is required
- Generally not configured for most routine operation



©2024 Waters Corporation


99

99

2424 Evaporative Light Scattering Detector

Waters™

- Detector gain
 - 1 ~ 1000
 - Will amplify peak and noise
 - Start with 200
- Sampling rate
 - 1 ~ 80Hz
 - HPLC method use 10 or 20 Hz
 - No significant improvement over 40 Hz
- Filter time constant
 - Normal is sufficient for general use
- Detector options:
 - Autozero is usually selected
 - Lamp on for acquisition, uncheck for shutdown method



©2024 Waters Corporation


100

100

2424 Evaporative Light Scattering Detector

Waters™

- Gas pressure
 - 20 ~ 60 PSI, or off
 - Have an effect on peak intensity
 - Generally starts with 40 PSI
- Nebuliser mode
 - Cooling, heating or off
 - Use only cooling
- Drift tube temperature
 - 5 to 100 °C
 - Higher temperature reduce baseline noise and peak intensity as this parameter is destructive
 - Alarm limit default at 25 °C is sufficient



©2024 Waters Corporation


101

101

2424 Evaporative Light Scattering Detector

Waters™

- Impact of drift tube temperature on the detection of dibutylphthalate
 - Above 45°C, DBP is nearly not detected
 - At 30°C, the baseline noise is 50% higher than at 45°C, but the response is 30 times greater
 - The best sensitivity is obtained with the lowest temperature value





©2024 Waters Corporation

102

102

2424 Evaporative Light Scattering Detector

Waters™



©2024 Waters Corporation


103

103

2424 Evaporative Light Scattering Detector

Waters™

- ELSD signal
 - Must be checked for recording signal
- ELSD signal - MBF
 - Similar to PDA MBF feature
 - Will smoothen the baseline drift so in the final result the baseline will appear as straight baseline
 - Will likely lost data of very small peaks



©2024 Waters Corporation


104

104

2424 Evaporative Light Scattering Detector

Waters™

- Events page is useful for time segmented gain adjustments
 - Some analytes concentration may be very high
 - Good to reduce gain to prevent oversaturating PMT
 - Likewise, gain can be increased to boost signal of less intense peak
 - Can also set Lamp to off
 - Can also set Gas Pressure to off
 - Both used for switching off ELSD after completing run



©2024 Waters Corporation


105

105

Sampling rate and Filter Time constants

Waters™

- General advice:
 - For UPLC application set filter time constant to “fast”
 - For sampling rate aim at 15-20 datapoints across the peak
- When sampling rate is set too high, the s/n will decrease:



©2024 Waters Corporation

106

106

General troubleshooting rule on LC hardware

Waters™

- First: Only change one thing at a time
- Second: Only call it a problem if it occurred twice
- Third: Put it back if it doesn't fix the issue
- Fourth: Throw bad parts away
- Fifth: Maintain a logbook for your LC system and write down important things
- Following these will reduce further downtimes – first time right
- John W. Dolan
 - <https://www.chromatographyonline.com/view/rules-thumb>
 - J.W. Dolan, LC•GC6(4), 304–308 (1988).

©2024 Waters Corporation

107

107

General LC operation consideration

Waters™

- Avoid adjusting pH of buffer after adding organics
 - Presence of organic affect pH
- Never top up buffer reservoir
 - Risk shifting pH and introducing debris into flowpath
- Use new column for LC method development
- If a column was used for ion-pairing methods, don't use on other methods
 - IP reagents will never ever be washed from column
- It is riskier to use column across multiple LC methods
- Not all C18 columns are equivalent
 - Use column coach to find out more about column selectivity
 - <https://find.waters.com/ColumnCoach/about>
- <https://www.chromatographyonline.com/view/seven-things-avoid-liquid-chromatography-laboratory-1>
- J.W. Dolan, LCGC North Am. 33(1), 18–22 (2015).

©2024 Waters Corporation

108

108

General LC operation consideration

Waters™

- Equilibrate and condition column properly
- Run blanks to prove system is clean before running samples
- It may be beneficial to run few prime injections

- <https://www.chromatographyonline.com/view/seven-things-avoid-liquid-chromatography-laboratory-1>
 - J.W. Dolan, LCGC North Am. 33(1), 18–22 (2015).
- <https://www.chromatographyonline.com/view/baker-s-dozen>
 - J.W. Dolan, LCGC North Am. 23(9), 1008–1012 (2005).

©2024 Waters Corporation

109

109

General LC column usage consideration

Waters™

- One column, one method
- In-line filters protect column from clogging effectively
- Clean column during idling
- Use temperature control – even for room temperature methods
- All column have void volume, make analytes elute after that
- Inject clean samples
- Use highest quality mobile phase available
- When validating methods, think of alternatives
- <https://www.chromatographyonline.com/view/lc-columns-top-10-list-1>
 - J.W. Dolan, LCGC North Am.21(3), 262–266 (2003).


©2024 Waters Corporation

110

110

Using Waters Knowledge Base

- Support.waters.com
- No need to signup/login
- If you have WKB (we call it KCS internally) article number, enter it in search box without WKB
- Or search with keywords
- Will return solution to general issues not available through Empower/MassLynx help files
- Contact us if your issue cannot be found on



©2024 Waters Corporation


111

111

Using Waters Knowledge Base

- Support.waters.com
- Searching with keywords works as well
- WKB articles cover informatics, instruments and consumables

Waters™



©2024 Waters Corporation

112

112

Waters™

Sample preparation – A very brief introduction

113

Why Do Sample Preparation?

Waters™

- Transfer analytes from matrix into a solution to be analyzed
- Remove matrix interferences
 - Increase signal to noise
 - Reduce variability in analytical results
 - Increase column lifetime
 - Reduces system downtime
 - Reduce chromatogram complexity
- Concentrate analyte
 - To meet the low detection limit
 - To stay within the sensitivity range of the available instrument

114

Sample preparation consideration

Waters™

- Ideally, HPLC/UPLC sample injections should be:
 - Clean, without visible solid particulates
 - Physical particulates will clog system
 - Optimal concentration to not contaminate instrument
 - Must account for matrix composition
 - Sample solvent must be compatible with mobile phase and system
 - Injecting neat strong solvent is generally advised against in Reversed Phase systems
 - Injection volume optimal for system and method
 - Not exceeding injection volume or consider the impact of injecting very large volumes

©2024 Waters Corporation


115

115

Common ways to prepare samples

Waters™

- Two major sample categories
 - Solid
 - Liquid
- At minimum (solids):
 - Grind
 - Weigh
 - Extract with solvent
 - Sonication
 - Filter through 0.22µm, OR, centrifuge/spin down
- <https://www.chromatographyonline.com/view/overview-sample-preparation>




©2024 Waters Corporation

116

116

Sample Preparation Techniques: Which One?

Waters™



©2024 Waters Corporation

117

117

Better ways to produce cleaner samples

Waters™

- Although you will get cleaner samples, these methods are:
 - More time consuming
 - Requires method development
 - Requires more consumables and equipment
 - Provides batch processing capability

©2024 Waters Corporation

118

118

Waters Oasis HLB SPE Products

Waters™

- Comes in both cartridges and well plate elution format
- Suitable for automated equipment or manual workflows



©2024 Waters Corporation


119

119

Some examples

Waters™

- General SPE will remove coloured contaminants and other matrix interferences
- Oasis PRiME effectively removed most coloured contaminants



©2024 Waters Corporation

120


120



Sample dilution consideration for LC analysis

Waters™

- Best to dilute in initial composition of mobile phase
- Avoid injecting in 100% organic solvent
 - Usually result in peak shape issue due to strong solvent effect
 - Right figure shows typical effect
- Sample concentration of around 1mg/mL is a good starting point, adjust when necessary to achieve good signal intensity
- For liquid samples, consider 100 or 1000 times dilution first



©2024 Waters Corporation

121

121

Waters™

Method transfer and development – A very brief introduction

122

Method transfer - Introduction

Waters™

- Why method transfer?
 - Adapting from official methods (ChP, GB, USP etc)
 - Reproducing method/results from other entity or published works
 - Obsolete column
 - Updated requirements, including accuracy and speed
- System where methods originate from are usually different from what you have
- System = LC + Column + Solvent + sample

©2024 Waters Corporation

123

123

Method transfer – criteria to consider

Waters™

- Validated method with example chromatograms
- Properties of donor and target system
 - Dwell volume and Mixer differences
 - Temperature control module, e.g. forced air or convection
 - Gradient shapes, i.e. if gradient curve was used
 - Detector settings, e.g. data rate, wavelength, bandwidth etc
- Preparation of mobile phase, online or offline mixing
- Steepness of gradient change, especially in UPLC methods
- Column chemistry, not all C18 are the same
- If transferring method from very old papers, consider modernisation

©2024 Waters Corporation

124

124

Required Information: Original Method, Results & Criteria

Waters™

- Column
 - Chemistry (ligand, particle size, brand)
 - Dimensions
- Conditions
 - Mobile phase
 - Flow rate
 - Gradient profile, including re-equilibration
 - Column temperature
- Sample
 - Diluent
 - Concentration
 - Injection volume
 - Molecular weight
- Chromatogram
 - Number of analytes
 - Retention required
 - Resolution required
- Quantitation
 - Limit of detection
 - Limit of quantitation
 - Linear dynamic range
 - Accuracy
 - Precision

©2024 Waters Corporation

125

125

Method transfer – desirable outcomes

Waters™

- Desirable outcomes from successful transfer:
 - No “missing” peaks
 - Equal or greater resolution between all critical pairs
 - Equal or greater signal to noise ratios
 - Equal or better precision for retention times, peak areas, peak heights and peak “shape” metrics (tailing, asymmetry, width)
 - Equal or greater dynamic range – upper and lower limits of quantitation meet or exceed current requirements
 - Order of elution is the same
 - Run time and cycle time are equal or faster
 - Solvent consumption is equal or less


©2024 Waters Corporation

126

126

Method transfer - Sources of System differences

Waters™



©2024 Waters Corporation

127

127

Empower 3 column calculator

Waters™

Columns Calculator

From... Describe your original method.

Column

Diameter (D):	4.600 mm
Length (L):	150 mm
Particle Size (dp):	5.0 µm
L/dp:	30,000

System

Dwell volume:	0.000 mL
---------------	----------

Method

Injection volume:	100 µL
Temperature:	30 °C
Run time:	10.00 min

To... Describe your target method.

Column

Diameter (D):	2.100 mm
Length (L):	50 mm
Particle Size (dp):	1.7 µm
L/dp:	29,412

System

Dwell volume:	0.000 mL
High pressure limit:	15,000 psi

Method

Flow rate:	Scaled: (0.613 mL/min)
Custom:	1.000 mL/min

Results:

Time (min)	Flow Rate (mL/min)	%A Water	%B Acetonitrile	%C Methanol	%D Water	Column Volumes
1.00	0.613	100.0	0.0	0.0	0.0	0.00

776 psi Maximum pressure

6,578 psi Maximum pressure

0.7 µL Injection volume

1.13 min Run time

©2024 Waters Corporation

128

128

Empower 3 column calculator

Waters™

- Column calculator comes with Empower 3
- Can be used to scale down or up a method
- Will take into account dwell volume of different systems
- Preferably using column of same chemistry in both systems
 - E.g. transferring Xbridge C18 from HPLC into Acquity BEH C18 in UPLC, same selectivity
- Use Waters Column Coach for looking for Waters column that provides comparable selectivity to your current column
 - <https://find.waters.com/ColumnCoach/about>

©2024 Waters Corporation

129

129


Method development

Waters™

- Why method development?
 - New product
 - Need a higher throughput LC method
 - Obsolete of column or previous LC system
 - End of life of LC method
- Core issue of LC method: resolution of all peaks in a chromatogram
 - What is resolution?
 - How much is enough?
 - Which factors will affect it?

$$R = \frac{2 \times (t_{R2} - t_{R1})}{W_1 + W_2} \text{ 或 } R = \frac{2 \times (t_{R2} - t_{R1})}{1.70 \times (W_{1,h/2} + W_{2,h/2})}$$

式中 t_{R2} 为相邻两色谱峰中后一峰的保留时间；
 t_{R1} 为相邻两色谱峰中前一峰的保留时间；
 W_1 、 W_2 及 $W_{1,h/2}$ 、 $W_{2,h/2}$ 分别为此相邻两色谱峰的峰宽及半高峰宽，见下图。



©2024 Waters Corporation

130


130

Resolution equation

Waters™

- k_2 : Retentivity
 - Chemical factor
 - Different mobile phase (MP), pH, composition of MP, temperature
- α : Separation factor
 - Chemical factor
 - Different stationary phase will improve resolution, or even change elution order
 - Derivatisation
- N : Efficiency
 - Theoretical plates
 - Physical/mechanical factor
 - Related to size of column and particles

$$R_s = \left(\frac{\sqrt{N_2^*}}{4} \right) \left(\frac{\alpha^* - 1}{\alpha^*} \right) \left(\frac{k_2^*}{1 + k_2^*} \right)$$



©2024 Waters Corporation

131

131

Method development – Other criteria

Waters™

- Sensitivity
- On-column loading
- Speed/throughput requirement
- Solvent consumption
- Ease of usage
- Column robustness/lifetime

©2024 Waters Corporation

132

132

Method development strategy

Waters™

- Look for reference, e.g. published papers or official methods
 - Must be soluble in aqueous or general reversed phase solvents
 - If there are no references usable, try scouting method
- General steps:
 - Scouting with 50 or 100 mm column
 - Run with standard first
 - Water/Acetonitrile is a good choice, adjust pH or use modifiers if required
 - Adjust k' for best retention
 - Adjust α for optimal selectivity
 - Consider different column length for best efficiency

©2024 Waters Corporation

133

133

Beneficial information during method development

Waters™

- Analyte solubility
- Number of analytes/peaks to be chromatographed
- Chemistry of analyte
- Type of detection required
- The effect of pH
- Detection level
- Matrix effect

©2024 Waters Corporation

134

134

Method development – choice of particles

Waters™

METHOD DEVELOPMENT FLEXIBILITY
WITH FOUR PREMIUM PARTICLE TECHNOLOGIES

BEH	CSH™	HSS	Solid-Core
BEH	CSH™	HSS	Solid-Core
- Good universal column choice for a wide variety of compounds - Exceptional peak shape for basic analytes at elevated pH - Stable across wide pH and temperature ranges	- Superior peak shape for basic compounds under acidic, low ionic strength conditions - Excellent MS performance with formic acid as a mobile phase modifier - Fast pH switching and column equilibration	- Increased retentivity over hybrid materials - Widest selectivity space with C₁₈, T3, C₁₈ SB, Cyano, and PFP chemistries - High strength silica (HSS) for mechanical stability	- Maximum efficiency - Increased sample throughput* - Lower column backpressure* * Compared to fully-porous particles of equivalent size.

©2024 Waters Corporation

135

135

Method development – choice of chemistry

Waters™

22 COLUMN CHEMISTRIES*

*For information on our line of bioseparation columns, visit: waters.com/biosep

BEH C₁₈	HSS T3	CORTECS C₁₈
BEH C₁₈ AX	HSS C₁₈	CORTECS C₁₈+
BEH C₈	HSS C₁₈ SB	CORTECS C₈
BEH Shield RP18	HSS PFP	CORTECS Phenyl Hexyl
BEH Phenyl	HSS Cyano	CORTECS T3
BEH HILIC		CORTECS Shield RP18
BEH Amide		CORTECS HILIC
BEH Z-HILIC		
CSH C₁₈		
CSH Fluoro-Phenyl		
CSH Phenyl-Hexyl	©2023 Waters Corporation. Waters, The Science of What's Possible, UPLC, XBridge, ACQUITY, CORTECS, and CSH are trademarks of Waters Corporation. All other trademarks are the property of their respective owners. January 2023 720006334EN KK-SIG	


©2024 Waters Corporation

136

136

Method development- choosing the optimal column length

Waters™



- L/d_p ratio is particularly useful when trying to determine which particle size packing material and column length may be necessary for a given application
- Easy – Simple assay, < 5 components
- Moderately challenging <10 components
- Difficult, general impurity profiling
- Extremely difficult, metabolites identification, toxicology screening


©2024 Waters Corporation

137

137

Method development- choosing the optimal column length

Waters™




- Key take-away:
 - Moving into UHPLC territory makes a great difference
 - Further cost savings and exceptional performance possible with UPLC
 - The key decision is throughput and economical consideration

©2024 Waters Corporation

138

138

Waters™



- Peak capacity (k') is a better measurement of gradient method performance
- Definition:
 - Theoretical number of peaks that can be separated in a given gradient time
 - See diagram to right for comparison between

©2024 Waters Corporation


139

139

Method transfer for efficiency gain

Waters™


How do you increase laboratory efficiency?



©2024 Waters Corporation

140

140



141